

---

**SOP NUMBER: 007**  
**TITLE: Thawing and Plating Cells**

---

Objective:

To thaw frozen cells so that they can be established in culture.

I. Materials:

A. Equipment

Biosafety cabinet  
Tabletop Centrifuge  
Freezer Liquid nitrogen  
Hemocytometer with cover glass  
Incubator, dual chamber  
Microscope, inverted  
Pipet-Aid  
Water bath @ 37°C

B. Supplies

Aspirating pipets: 2 ml  
Cryogloves  
Dry ice  
Face shield  
Nitrile gloves  
Plastic Serological Pipettes: 1ml, 2 ml, 5,10 ml  
Polystyrene centrifuge tubes with screw cap, 15 ml

C. Media and Chemicals

1. Chemicals  
Growth medium of cell line to be thawed  
Phosphate Buffered Saline w/o calcium and magnesium (1x)  
Trypan blue stain 0.4%

II. Method:

1. Label 15 ml conical tube with cell line designation. Pipet 10 mls of PBS into each tube.
3. Remove frozen vial of cells from the liquid nitrogen freezer and place on dry ice.
4. Thaw cells rapidly in a 37°C water bath while swirling. The swirling prevents formation of intracellular ice. **Safety Note:**The operator should wear safety glasses and a face shield.
5. After thawing, slowly add 1 ml of PBS drop-wise to the vial of cells. This is to reduce osmotic swelling upon rehydration and improve post-freeze viability.
6. Pipet the 2 mls of cell suspension from the cryovial and place in the 15 ml conical tube containing PBS.
7. Wash and pellet cells at 150 Xg to remove cryopreservative.
8. Count cells and determine post-freeze viability with Trypan Blue dye. See SOP #290: Cell Count Determination. If post-freeze viability is low, cells should be reseeded at a high concentration to increase chances of recovery.

---

**SOP NUMBER: 007**  
**TITLE: Thawing and Plating Cells**

---

9. Establish growth in pre-warmed culture medium that has been tested for sterility.

Reference: Margaret B. Penno, "Cryopreservation" pp108-112, Basic Cell Culture-A Practical Approach, 1998.

---

The Johns Hopkins Medical Institutions  
The Genetic Resources Core Facility: The Cell Center  
Blalock 1016/600 N. Wolfe Street Baltimore, Maryland 21287  
Director: Margaret B. Penno, Ph.D. Tel: 410-614-0060 E-mail: penno@jhmi.edu

---